





SMALLER DEEPER BETTER

NanoFx for enhanced MicroFx



STEP 1
Assemble and orient on lesion



STEP 2
Gently tap to depth stop



STEP 3:
Push Thumble to extract
PluriStik. Repeat.



*To order Tel:
01793 575050
see product
details and ordering
information below*

CAT NO:	DESCRIPTION
500-1020	NanoFx Handle - Knee
FURS-2101	5 PleuriStik™ - Knee
ROCK-PL1	1 PleuriStik™ - Knee
FURS-0100	5 Thumbles

CAT NO:	DESCRIPTION
5500-4010	NanoFx Handle - Talus
FURS-4101	5 PleuriStik™ - Talus
ROCK-AC1	1 PleuriStik™ - Talus
ROCK-TM1	1 Thumble

JOINT OPERATIONS

Suite 20, Anchor Business Centre, Frankland Road, Blagrove, Swindon, Wiltshire SN5 8YZ

T 01793 575050 F 01793 575040 E service@jointoperations.co.uk

