



arthrosurface®

ARTICULATION. FIXATION. PRESERVATION.

The Arthrosurface® ToeMotion™ System restores mobility and maintains native biomechanics using a dual curved HemiCAP DF® and a new modular tray-style phalangeal implant with a threaded baseplate. Fourth generation fixation components provide stable constructs on both sides of the joint.



www.arthrosurface.com



- Inlay arthroplasty provides stability and preserves bone
- Double rock-solid fixation with morse taper interlock
- Long-term clinical history with over 20,000 MTP implants
- Single tray system with multiple implant options