HealFix Ti Suture Anchors Family

Various options in anchor design, diameter, suture size/color, with/without needles

reddot winner 2020

The Red Dot Award:Product Design only awards this sought-after seal of quality to products that feature an outstanding design. In 2020, HealFix Suture Anchor won the Red Dot award in the category"Medical Devices and Technology".



1.3mm with needles
2.0mm with needles
2.4mm with needles



2.0mm with needles
2.0mm



2.8mm with needles
2.8mm



3.5mm with needles
3.5mm



5.0mm with needles
5.0mm
6.5mm with needles
6.5mm



4.5mm Plus
5.5mm Plus
6.5mm Plus





rejoin

Hangzhou Rejoin Mastin Medical Device Co., Ltd

Headquarter Office: 22 Xinyan Rd, Linping Economic Development Zone, Hangzhou, Zhejiang, 311100, P.R. of China
Beijing Branch: Room 1605, Jinghui building, No. 118 Jianguo Road, Chaoyang District, Beijing, P.R. of China
Shanghai Branch: 309A, L1-B, Hongqiao World Center, 1588 Zhuguang Rd, Qingpu District, Shaghai, P.R. of China

contact@rejoin.com **@rejoinmedical** **www.rejoin.com**
www.rejoin-medical.com

The content in this document is provided for reference only

rejoin


HealFix Titanium Anchor Fixation System

Wide Range of Applications

More options for sizes, threads & sutures

Pre-loaded with WilSuture Ultra sutures


Ergonomic handle design provides surgeons comfort and ease of use



HealFix Ti Suture Anchor

- High-Low double thread design enhances the pull-out strength
- Two independent eyelets provide optimal fixation of both sutures & eliminate suture twining
- Sharpened tip facilitates easy insertion





HealFix Plus Ti Suture Anchor

- Fully Helix transitional thread design increases contact areas with cortical bone for high pull-out strength
- Built-in suture hang-up hole eliminates Cheese-Wire effect between sutures and bone after screw implantation
- Self-punching design makes implantation easier



HealFix Mini Ti Suture Anchor

- NiTi arcs design provides subcortical purchase with outstanding pull-out strength
- Smaller size allows comprehensive fixation and restoration of natural footprint
- 1.3mm minimum size provides optimized soft tissue fixation for DEX applications

Ordering Information

HealFix

Code	Description
800500600	HealFix Ti Suture Anchor, 2.0mm, w/two #2-0 sutures W&WB
800500800	HealFix Ti Suture Anchor, 2.0mm, w/two #2-0 sutures W&WB, w/needles
800501400	HealFix Ti Suture Anchor, 2.8mm, w/#2 suture WB
800501200	HealFix Ti Suture Anchor, 2.8mm, w/two #2-0 sutures W&WB
800501000	HealFix Ti Suture Anchor, 2.8mm, w/two #2-0 sutures W&WB, w/needles
800503500	HealFix Ti Suture Anchor, 3.5mm, w/two #2 sutures W&WB
800503300	HealFix Ti Suture Anchor, 3.5mm, w/two #2 sutures W&WB, w/needles
800503900	HealFix Ti Suture Anchor, 5.0mm, w/two #2 sutures W&WB
800503700	HealFix Ti Suture Anchor, 5.0mm, w/two #2 sutures W&WB, w/needles
801310000	1.6mm Drill for 2.0mm HealFix
801310100	Drill Gudie for 2.0mm HealFix
801310200	1.8mm Drill for 2.8mm HealFix
801310300	Drill Guide for 2.8mm/3.5mm HealFix
801310400	2.5mm Drill for 3.5mm/5.0mm HealFix
801310600	3.5mm Drill for 6.5mm HealFix

HealFix Plus

Code	Description
800504500	HealFix Plus Ti Suture Anchor, 4.5mm, w/two #2 sutures W&WB
800504700	HealFix Plus Ti Suture Anchor, 5.5mm, w/two #2 sutures W&WB
801310800	3.5mm Drill for 4.5mm HealFix Plus
801310900	4.5mm Drill for 5.5mm HealFix Plus
801311000	5.5mm Drill for 6.5mm HealFix Plus

HealFix Mini

Code	Description
800505300	HealFix Mini Ti Suture Anchor, 1.3mm, w/3-0# suture W, w/tapercut needles
800505400	HealFix Mini Ti Suture Anchor, 1.3mm, w/3-0# suture WB, w/tapercut needles
800505900	HealFix Mini Ti Suture Anchor, 1.8mm, w/2-0# suture W, w/tapercut needles
800506000	HealFix Mini Ti Suture Anchor, 1.8mm, w/2-0# suture WB, w/tapercut needles
801911100	Drill for 1.3mm HealFix Mini
801911200	2.0mm Drill for 1.8mm HealFix Mini